Request license update

Åbn CodeMeter Control Center



Jeres Dongle ID vil stå i boksen til venstre. Klik på "License Update"

ScodeMeter Control Center			×
File Process View Help			
License Events Borrow			
GOM License 2-2951802	Name:	GOM License	ø
	Serial:	2-2951802	9
	Version:	CmStick 2.02	ଷ
	Capacity:	90 % free (352512 Bytes)	
	Status:	C S Disabled C S Disabled	
	renne lindata	Change Descured	
CodeMeter is started.	Lerise opdate	cjett Unange Password	WebAdmin

Klik "Next"





Vælg "Create license request" – "Next"



Vælg "Extend existing license" - "Next"



Der skulle gerne være flueben ved GOM Klik "Next"



Vælg hvor filen skal gemmes – klik "Commit"



Klik "Finish"



Nu ligger der en WIBU Control fil, der hvor du har specificeret.



Upload filen til GOM via GOM Dongle Update Form https://www.gom.com/services/licence-update.html

License Update

GOM Software 2020 requires a new license. To update your license dongle to the new GOM Software 2020, please use the Dongle Update Form.

Dongle Update to GOM Software 2020			
Name*:	Please insert your name		
Email*:	Please insert your email		
Email repeat*:	Please insert your email again		
Comment:	Add your comment here		
Dongle License Request File*:	Vælg fil Der er ikke valgt nogen fil		
Antispam protection*:	Please insert the code below		
	v562mW		
	reset submit		
* required			

Herefter modtager du en opdateret licensfil fra GOM, som skal læses ind på jeres dongle.

Update dongle license

Åbn CodeMeter Control Center



Jeres Dongle ID vil stå i boksen til venstre. Klik på "License Update"

S CodeMeter Control Center	×
File Process View Help	
License Events Borrow	
GOM License Name: GOM License 2-2951802	Ø
Serial: 2-2951802	9
Version: CmStick 2.02	ଷ
Capacity: 90 % free (352512 Bytes)	
Status: 🔿 😋 Disabled	
Generation Generation	
License Update Eject Change Password	
CodeMeter is started.	WebAdmin

Klik "Next"





Vælg "Import license update" – "Next"



Vælg den opdaterede licensfil fra GOM – klik "Commit"

